

AESTHETIC MEDICINE



GUNA®-COLLAGEN is a sterile solution made with isotonic sodium chloride solution.

The strength at 6X supports the collagen physiological function, as the homeopathic biological substrate acts as supporter of the endogenous production of collagen. The local mesotherapy technique of administration targets the area of activity.

GUNA®-COLLAGEN shows a “flattening” effect on wrinkles, that improves as the protocol goes on; the effect is due to a better metabolism on the area of application not to a filling activity.

USES

Skin tone improvement, treatment of face and neck wrinkles, stimulation of connective tissue metabolism.

INGREDIENTS

Active ingredients:

Collagen, Porcine 6X.

Inactive ingredient:

Sterile Isotonic Sodium Chloride Solution.

DOSAGE FORMS

2 ml glass ampoules in package of 10 ampoules
2 ml glass ampoules in package of 50 ampoules

DIRECTIONS

Standard protocol for i.m. administration:

Skin Tone Improvement

1 ampoule 1-3 times a week according to severity and clinic evolution

Standard protocol according to mesotherapy technique:

Treatment of Face and Neck Wrinkles

1 ampoule / face . Using a 4 mm, 30G needle, make the classic intradermal injection according to the mesotherapy technique.

Stimulation of connective tissue metabolism:

Treatment of Neck Wrinkles with Mesotherapy Technique:

0.3 ml. Using a 13 mm, 30G needles, cannulate the wrinkle by moving the needle gently to left and right, while injecting the content of the syringe.

Mesotherapy Technique Standard Protocol consists of one session a week for 5-7 weeks. For prolonged treatments, 2 sessions for the first week, 1 session a week for 1 month and then 1 session a month are recommended.